

Hg analysis for solid samples

Oulu laboratory


Direct mercury analysis for solid matrix (fish, plants, sludge).

NIC MA-3000 analyser with direct combustion and cold vapour atomic absorption spectrometry (CV-AAS).


- Analysis based on US-EPA 7473 standard
- Accredited for fish matrix
- Limit of detection for fish 0.02 mg/kg
- Fish preparation service also available
- Contact persons: Mika Sarkkinen and Ritva Nilivaara (firstname.lastname@syke.fi)



© Katri Salonen



© Katri Salonen



© Mika Sarkkinen