

Centrifuge for harvesting microalgae

gear gallery

Centrifuge for separating microalgae or other particles from suspension.

- AlfaLaval Clara 20 separator for application which require high hygienic conditions
- Can process 100-500 liters of algal culture in an hour
- In addition to the high separation efficiency, gentle treatment of the biomass,
- Robust and easy to operate
- For technical details, see the product <u>leaflet</u> on <u>manufacturer's web site</u>.
- Contact Kristian Spilling (firstname.surname@syke.fi)



Finnish Environment Institute

